

# Semi-automated Protein Precipitation Method Using Tip-on-Tip Technology and the Integra ASSIST PLUS for Analysis of Benzodiazepines in Biological Fluids

AAFS 2020 Poster Summary: Dani Mata, Orange County Crime Lab, Santa Ana, CA

**HIGHLIGHTS:** Replaces manual vortex and centrifugation



## INTRODUCTION AND OBJECTIVE

- Introduce innovative, automated INTip™ solution as an alternative to centrifugation utilizing a patent-pending Tip-on-Tip (ToT) technology
- Low Porosity Filtration Tips and an Integra ASSIST PLUS minimize laborious manual steps such as vortex mixing and centrifugation
- Extend the use of the Integra for the extraction method that Orange County Crime Lab (OCCL) uses to analyze benzodiazepines. The extraction enables the OCCL to quantitate 23 benzodiazepines, z-drugs and antihistamines, and qualitatively detect 7 metabolites in urine, blood and tissue specimens.



Table 2. Tip-on-Tip protein precipitation procedure and sample cleanup using the Integra ASSIST PLUS

## METHODS

Table 1. Current manual protein precipitation procedure using vortex and centrifugation

<b>Protein Precipitation</b>	<ul style="list-style-type: none"> <li>• Manually pipette 250 µL sample, 50 µL internal standards, 750 µL acetonitrile (ACN)</li> <li>• Vortex and Centrifuge</li> <li>• Manually decant supernatant to new test tube</li> </ul>
<b>Sample Cleanup</b>	<ul style="list-style-type: none"> <li>• Aspirate/dispense supernatant 3x using WAX-XTR tips</li> </ul>
<b>Prep for injection</b>	<ul style="list-style-type: none"> <li>• Manually transfer 50 µL of eluant, add 800 µL initial mobile phase into LC vials</li> </ul>

<b>Protein Precipitation</b>	<ul style="list-style-type: none"> <li>• Manually pipette 250 µL sample/standard in test tube. Add 50 µL of internal standard</li> <li>• Add 750 µL of cold ACN</li> <li>• Mix sample and ACN using wide bore tips - crash proteins</li> <li>• Attach wide bore tips to Low Porosity Filtration tips and dispense supernatant into new test tubes</li> </ul>
<b>Sample Cleanup</b>	<ul style="list-style-type: none"> <li>• Aspirate/dispense supernatant 3x using WAX-XTR tips</li> </ul>
<b>Prep for injection</b>	<ul style="list-style-type: none"> <li>• Manually transfer 50 µL of eluant, add 800 µL initial mobile phase into LC vials</li> </ul>

# WORKFLOW USING INTEGRA ASSIST PLUS




Figure 1. Schematic represents workflow for ToT Filtration used for protein precipitation and WAX-XTR cleanup method.

## RESULTS

- **Casework** - Percent differences of DPX Tip-on-Tip method versus manual method for AM/PM blood, urine, liver, brain and gastric casework ranged from -15.5% to 15.4%.
- **Reproducibility** - Concentrated liver homogenate had within-run precision of 2.5% and 7.9% for alprazolam and hydroxyzine, respectively.
- **Linearity** -  $R^2$  was greater than 0.99 for all drugs.
- **Accuracy** - As shown in figure 2, QC's of porcine blood were extracted with the manual and automated DPX extraction methods. Percent differences ranged from -5 to 13%




Figure 2. Percent Difference for the High and Low QC between the Manual Extraction and Automated Extraction. QC's are made in porcine blood.

## CONCLUSIONS

These initial investigative results demonstrate that protein precipitation can be automated for a SWGTOX validated forensic application of blood, urine and tissues by LC-MS/MS. Adaptation of an automated crash method for forensic laboratories offers several advantages including reducing the hands-on time requirement for the analyst, minimizing potential human errors, and further negating the long-term health effects of repeated pipetting.