

Innovative NiXTip™ Technology for Efficient Size Selection in Sequencing Library Preparation

SLAS 2025 Poster Summary: Paul F. Meeh, Matt G. Fitts, and William E. Brewer

HIGHLIGHTS: Fully automated, fast (< 30 min), customizable



INTRODUCTION

Solid Phase Reversible Immobilization (SPRI) with magnetic beads has been the gold standard for DNA clean up and size selection in Next Generation Sequencing (NGS) protocols. The use of magnetic beads eliminates the need for centrifugation but has its own inherent draw backs such as bead capture time, bead carryover, bead resuspension and the requirement for magnets. DPX Technologies has developed NiXTips™ that eliminate magnetic beads and magnets for DNA clean up and size selection in NGS library preparations. NiXTips incorporate a proprietary silica media within a pipette tip to utilize aspirate and dispense steps, which streamline DNA recovery while maximizing efficiency. Applications include DNA capture from PCR and enzymatic reaction cleanup, as well as double-sided size selection for sequencing libraries. Results demonstrate almost 100% recovery of DNA depending on target fragment size and buffer composition of the sample. The data presented herein shows that NiXTips simplify workflows while providing high yield of DNA targets.

MATERIALS AND METHODS



Figure 1. NiXTips shown on a Hamilton® Robotics platform.

NiXTips were utilized with the accompanying proprietary room-temperature binding buffer for optimal performance (DPX Technologies, Columbia, SC). Wash buffer (80% ethanol) and elution buffer (molecular biology grade water) were also utilized but are not provided in the NiXTips kit. DNA samples prepared using an in-house PCR product with an amplicon size of 205 bp were used for the cleanup and Sanger sequencing experiments. DNA molecular weight marker (MWM) VIII with a size range of 19 to 1114 bp (Millipore Sigma, Burlington, MA) was used for the size selection and binding cycle studies. PCR reactions were performed on a QuantStudio™ 3 (Thermo Fisher Scientific, Waltham,

MA) with the TaqMan™ Genotyping MM from Thermo Fisher Scientific. For these studies, 25 µL of DNA sample was used. All samples, reagents, and tips were loaded onto a Hamilton® NIMBUS automated liquid handler (ALH). Multiple iterations were performed with various binding buffer ratios and binding cycles to perform size selection and optimize yield, respectively.

1	Sample Prep	ALH picks up transfer tips and aspirates binding buffer to achieve the desired binding buffer to sample ratio.
2	Bind	Binding buffer is added to the sample and mixed by aspirating and dispensing 20 times.
3	Wash	NiXTips aspirate and dispense sample up to 40 times to bind DNA.
4	Elute	Wash with 80% ethanol by aspirating and dispensing 50 µL of wash buffer with the NiXTip, 1 time.
		Dry the proprietary silica media using an air blowout to eliminate any residual ethanol.
		NiXTips aspirate and dispense 27 µL of elution buffer at 56°C, 2 times.

Table 1. NiXTip Method. For elution, 27 µL is used to achieve a final eluent volume of 25 µL due to a 2 µL loss to the media.

The eluted DNA was analyzed to assess purity and yield using a Thermo Fisher Scientific NanoDrop One Microvolume UV-Vis Spectrophotometer and an Agilent 2100 Bioanalyzer (Santa Clara, CA). To confirm size selection, agarose gel electrophoresis was performed. Compatibility for downstream applications was verified by Sanger sequencing of the GAPDH gene.




Figure 2. Schematic of NiXTips workflow.

