




Isolation of gDNA from Cell Lines/Cultures



Technology: Proprietary Silica Media
Kits include: Tip Format of Choice, Proprietary Lysis, Wash (1 & 2), and Elution Buffers
Applications: Adherent or suspension cell lines for downstream PCR, RT-PCR or ddPCR

INTEGRA



125 µL

Capacity	# of Cells	Sample Volume	Elution Volume	Processing Time
14 µg	<100,000	25-100 µL	25-100 µL	18 min

HAMILTON



300 µL

Capacity	# of Cells	Sample Volume	Elution Volume	Processing Time
11 µg	<100,000	25-100 µL	25-100 µL	30 min



1 mL

Capacity	# of Cells	Sample Volume	Elution Volume	Processing Time
43 µg	<2,000,000	25-250	50-250 µL	36 min

Coming Soon: 200 µL Tecan, 300 µL Integra, 1250 µL Integra